

SMARTX: Porewater Ammonium

May 2026 Samples

2026-09-03

Contents

0.1	Assessing standard Curves	4
0.2	Update Standard Log	4
0.3	Dilution Corrections - ensure the latest dilution is kept	6
0.4	Performance Check	6
0.5	Analyze the Check Standards	7
0.6	Analyze Blanks	8
0.7	Analyze Duplicates	9
0.8	Spikes	10
0.9	Matrix Effects	11
0.10	Create Log for QAQC	12
0.11	Merge samples with Metadata & check all samples are present	15
0.12	Visualize Data	16
0.13	Organize	20

##Run Information

```
cat("Run Information: ") #lets you know what section you're in
```

Run Information:

```
library(here)
```

here() starts at C:/Users/blkleybd/OneDrive - Smithsonian Institution/Documents/GitHub/GCREW

```
#set the run date & user name
```

```
sample_year <- "2026"
```

```
sample_month <- "May"
```

```
user <- "Brenden Blakley"
```

```
#identify the files you want to read in
```

```
#read in as a list to accommodate multiple runs in a month
```

```
NH3_files <- c("Raw Data/20260624_SMARTX_Run1_May2026.csv",  
              "Raw Data/20260626_SMARTX_Run2_May2026.csv",  
              "Raw Data/20260807_SMARTX_Run3_May2026.csv",  
              "Raw Data/20260626_2M_SMARTXrr_May2026.csv"  
              )
```

```
# Define the file path for QAQC log file - NO Need to change just check year
```

```
file_path <- here("Porewater QAQC Logs", "GCREW_NH4_QAQC_Log.csv")
```

```
file_path_QAQC <- here("Porewater QAQC Logs", "GCREW_SEAL_Check_QAQC_Log.csv")
```

```
final_path <- "Processed Data/GCREW_SMARTX_Porewater_NH4_202605.csv"
```

```
#record any notes about the run or anything other info here:
```

```
run_notes <- "SMARTX_441-80 Not Collected; One PE Check Passes One Fails on 20260626, likely air bubble"
```

```
#Add SMARTX Prefix if necessary to samples
```

```
Add_Prefix <- "YES"
```

```
#Collection_Dates = "Raw Data/Sample_Collection_Dates_2024.csv"
```

```
cat(run_notes)
```

SMARTX_441-80 Not Collected; One PE Check Passes One Fails on 20260626, likely air bubble

##Setup


##Read in metadata and create similar sample IDs for matching to samples

```
##Import Data & Clean
```

0.1 Assessing standard Curves

```
## 'geom_smooth()' using formula = 'y ~ x'
```


NH3 Standard Curve




Test	Run_Date	R2_Flag	R2	Slope
Ammonia 2	06/24/2026	R2_Pass	0.9998859	1.893843
Ammonia 2	06/26/2026	R2_Pass	0.9998100	1.856283
Ammonia 2	08/07/2026	R2_Pass	0.9985244	1.934832

0.2 Update Standard Log

```
## 'summarise()' has grouped output by 'Test', 'Intercept', 'Slope', 'R_Squared',  
## 'User'. You can override using the '.groups' argument.
```





0.3 Dilution Corrections - ensure the latest dilution is kept

```
## [1] "Reran Samples: c(\"321-20\", \"323-20\", \"331-20\", \"332-20\", \"333-20\", \"633-20\", \"312-20\")"
```

```
## [1] "Diluted Samples: c(\"321-20\", \"323-20\", \"331-20\", \"332-20\", \"333-20\", \"633-20\", \"312-20\")"
```


```
## [1] "No Naming Issues Detected"
```

0.4 Performance Check

Test	Run_Date	PE_Flag	PE_Conc	PE_Target_Conc
Ammonia 2	06/24/2026	Performance Check Within 25% - PROCEED	1.086850	1.034
Ammonia 2	06/26/2026	Performance Check Within 25% - PROCEED	1.272801	1.034
Ammonia 2	06/26/2026	Performance Check Within 25% - PROCEED	1.272801	1.034
Ammonia 2	08/07/2026	Performance Check Within 25% - PROCEED	1.095402	1.034

0.5 Analyze the Check Standards


Test	Run_Date	RSV_Flag	RSV	RSV_Cutoff
Ammonia 2	06/24/2026	RSV WITHIN RANGE - PROCEED	0.0242693	0.25
Ammonia 2	06/26/2026	RSV WITHIN RANGE - PROCEED	0.0242693	0.25
Ammonia 2	08/07/2026	RSV WITHIN RANGE - PROCEED	0.0242693	0.25




Test	Run_Date	CHK_Flag
Ammonia 2	06/24/2026	>60% of Checks Pass - PROCEED
Ammonia 2	06/26/2026	>60% of Checks Pass - PROCEED
Ammonia 2	08/07/2026	>60% of Checks Pass - PROCEED

0.6 Analyze Blanks

Test	Run_Date	BLK_Pct_Flag	Mean_Blkc_Conc	Quantile_25
Ammonia 2	06/24/2026	>60% of Blanks Pass - PROCEED	-0.0001170	0.2514288
Ammonia 2	06/26/2026	>60% of Blanks Pass - PROCEED	-0.0015411	1.9243220
Ammonia 2	08/07/2026	>60% of Blanks Pass - PROCEED	-0.0484810	3.3251115




0.7 Analyze Duplicates



Test	Run_Date	Dup_Flags
Ammonia 2	06/24/2026	>60% of Dups Pass - PROCEED
Ammonia 2	06/26/2026	>60% of Dups Pass - PROCEED
Ammonia 2	08/07/2026	>60% of Dups Pass - PROCEED


0.8 Spikes



Test	Run_Date	Spike_Flags
Ammonia 2	06/24/2026	>60% of Spikes have a Diff <50% - PROCEED
Ammonia 2	06/26/2026	>60% of Spikes have a Diff <50% - PROCEED
Ammonia 2	08/07/2026	>60% of Spikes have a Diff <50% - PROCEED
NA	NA	>60% of Spikes have a Diff <50% - REASSESS

0.9 Matrix Effects


```
## Warning in left_join(Matrix, Spike_ready, by = join_by(Date_Num)): Detected an unexpected many-to-many
## i Row 7 of 'x' matches multiple rows in 'y'.
## i Row 166 of 'y' matches multiple rows in 'x'.
## i If a many-to-many relationship is expected, set 'relationship =
## "many-to-many"' to silence this warning.
```



Test	Run_Date	Matrix_Flags
Ammonia 2	06/24/2026	Matrix Has Diff <50% - PROCEED
Ammonia 2	06/26/2026	Matrix Has Diff <50% - PROCEED
Ammonia 2	08/07/2026	Matrix Has Diff <50% - PROCEED

##Unit Conversion

```
##Sample Flagging - Within range of standard curve
```




```
##Add QAQC Flags
```


0.10 Create Log for QAQC

```
## Warning in left_join(Top_Standards, PE_Check, by = join_by(Test, Run_Date)): Detected an unexpected
## i Row 2 of 'x' matches multiple rows in 'y'.
## i Row 2 of 'y' matches multiple rows in 'x'.
## i If a many-to-many relationship is expected, set 'relationship =
##   "many-to-many"' to silence this warning.


## 'summarise()' has grouped output by 'Test', 'Run_Date', 'Top_Std_Absorbance',
## 'Mean_PE_Conc'. You can override using the '.groups' argument.

## Warning: Removed 6 rows containing missing values or values outside the scale range
## ('geom_point()').
```





```
## Warning: Removed 1 row containing missing values or values outside the scale range
## ('geom_point()').
```



```
##Write QAQC Log
```

0.11 Merge samples with Metadata & check all samples are present

```
## Merge Samples with Metadata
```


```
## Some sample IDs are missing from data.
```

```
## [1] "SMARTX_441-80"
```


0.12 Visualize Data

```
## Visualize Data
```



C3: Porewater NH3




C4: Porewater NH3



C3: Porewater NH3 by Replicate





0.13 Organize

#Export Data

#end