

SMARTX Sulfide Microplate Template

SMARTX 2026 July

Name: Melanie

Plate ID: plate 1

Date: 7/20/26

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S ⁶¹⁰⁻ ₂₀	S	S	S ³¹⁰⁻ ₁₂₀	S	S	S ³¹¹⁻ ₁₂₀	S	S
B	Std 1	Std 1	Std 1	S ⁶¹⁰⁻ ₄₀	S	S	S0	S0	S0	S ³¹²⁻ ₁₀	S	S
C	Std 2	Std 2	Std 2	S ⁶¹⁰⁻ ₈₀	S	S	S ³¹⁰⁻ _{Dup 120}	Dup	Dup	S ³¹²⁻ ₂₀	S	S
D	Std 3	Std 3	Std 3	S ⁶¹⁰⁻ ₁₂₀	S	S	S3	S3	S3	S ⁶¹⁰⁻ _{Spike 10}	Spike	Spike
E	Std 4	Std 4	Std 4	S ³¹⁰⁻ ₁₀	S	S	S ³¹¹⁻ ₁₀	S	S	S ³¹¹⁻ ₁₀ Spike	Spike	Spike
F	Std 5	Std 5	Std 5	S ³¹⁰⁻ ₂₀	S	S	S ³¹¹⁻ ₂₀	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S ³¹⁰⁻ ₄₀	S	S	S ³¹¹⁻ ₄₀	S	S	S ³¹²⁻ _{Dup 20}	Dup	Dup
H	S ⁶¹⁰⁻ ₁₀	S	S	S ³¹⁰⁻ ₈₀	S	S	S ³¹¹⁻ ₈₀	S	S	S4	S4	S4

12:09

X100

SMARTX Sulfide Microplate Template

SMARTx 2026 July

Name: my name

Plate ID: plate 2

Date: 7/20/26

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S ³¹²⁻ ₈₀	S	S	S ⁶¹³⁻ ₂₀	S	S	S ⁶²⁰⁻ ₂₀	S	S
B	Std 1	Std 1	Std 1	S ³¹²⁻ ₁₂₀	S	S	S0	S0	S0	S ⁶²⁰⁻ ₄₀	S	S
C	Std 2	Std 2	Std 2	S ³¹³⁻ ₁₀	S	S	S ⁶¹³⁻ _{Dup 20}	Dup	Dup	S ⁶²⁰⁻ ₈₀	S	S
D	Std 3	Std 3	Std 3	S ³¹³⁻ ₂₀	S	S	S3	S3	S3	S ³¹²⁻ _{Spike 40}	Spike	Spike
E	Std 4	Std 4	Std 4	S ³¹³⁻ ₄₀	S	S	S ⁶¹³⁻ ₄₀	S	S	S ⁶¹³⁻ _{Spike 40}	Spike	Spike
F	Std 5	Std 5	Std 5	S ³¹³⁻ ₈₀	S	S	S ⁶¹³⁻ ₁₈₀	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S ³¹³⁻ ₁₂₀	S	S	S ⁶¹³⁻ ₁₂₀	S	S	S ⁶²⁰⁻ _{Dup 80}	Dup	Dup
H	S ³¹²⁻ ₄₀	S	S	S ⁶¹³⁻ ₁₀	S	S	S ⁶²⁰⁻ ₁₀	S	S	S4	S4	S4

13:22

X100

SMARTX Sulfide Microplate Template

Name: Melanie

Date: 7/20/26

SMARTX 2026 July
Plate ID: plate 3

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S 320- 40	S	S	S 321- 80	S	S	S 322- 80	S	S
B	Std 1	Std 1	Std 1	S 320- 20	S	S	S0	S0	S0	S 322- 120	S	S
C	Std 2	Std 2	Std 2	S 320- 40	S	S	Dup 321- 80	Dup	Dup	S 323- 10	S	S
D	Std 3	Std 3	Std 3	S 320- 80	S	S	S3	S3	S3	Spike 620- 120	Spike	Spike
E	Std 4	Std 4	Std 4	S 320- 120	S	S	S 321- 120	S	S	Spike 321- 120	Spike	Spike
F	Std 5	Std 5	Std 5	S 321- 10	S	S	S 322- 10	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S 321- 20	S	S	S 322- 20	S	S	Dup 323- 10	Dup	Dup
H	S 620- 120	S	S	S 321- 40	S	S	S 322- 40	S	S	S4	S4	S4

14:04

X100

SMARTX Sulfide Microplate Template

SMARTX 2026 July

Name: melanre

Plate ID: plate 4

Date: 7/20/26

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S ³²³⁻ ₄₀	S	S	S ⁶³⁰⁻ ₁₀	S	S	S ³³⁰⁻ ₁₀	S	S
B	Std 1	Std 1	Std 1	S ³²³⁻ ₈₀	S	S	S0	S0	S0	S ³³⁰⁻ ₂₀	S	S
C	Std 2	Std 2	Std 2	S ³²³⁻ ₁₂₀	S	S	Dup ⁶³⁰⁻ ₁₀	Dup	Dup	S ³³⁰⁻ ₄₀	S	S
D	Std 3	Std 3	Std 3	S ⁶²³⁻ ₁₀	S	S	S3	S3	S3	Spike ³²³⁻ ₂₀	Spike	Spike
E	Std 4	Std 4	Std 4	S ⁶²³⁻ ₂₀	S	S	S ⁶³⁰⁻ ₂₀	S	S	Spike ⁶³⁰⁻ ₂₀	Spike	Spike
F	Std 5	Std 5	Std 5	S ⁶²³⁻ ₄₀	S	S	S ⁶³⁰⁻ ₄₀	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S ⁶²³⁻ ₈₀	S	S	S ⁶³⁰⁻ ₈₀	S	S	S ³³⁰⁻ ₄₀ Dup	Dup	Dup
H	S ³²³⁻ ₂₀	S	S	S ⁶²³⁻ ₁₂₀	S	S	S ⁶³⁰⁻ ₁₂₀	S	S	S4	S4	S4

14:55

X100

SMARTX Sulfide Microplate Template

SMARTX 2026 July

Name: Melanie

Plate ID: Plate 5

Date: 7/20/26

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S 330-120	S	S	S 332-40	S	S	S 333-40	S	S
B	Std 1	Std 1	Std 1	S 331-10	S	S	S0	S0	S0	S 333-80	S	S
C	Std 2	Std 2	Std 2	S 331-20	S	S	S 332-40 Dup	Dup	Dup	S 333-120	S	S
D	Std 3	Std 3	Std 3	S 331-40	S	S	S3	S3	S3	Spike 330-80	Spike	Spike
E	Std 4	Std 4	Std 4	S 331-80	S	S	S 332-80	S	S	Spike 332-80	Spike	Spike
F	Std 5	Std 5	Std 5	S 331-120	S	S	S 332-120	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S 332-10	S	S	S 333-10	S	S	Dup 333-120	Dup	Dup
H	S 330-80	S	S	S 332-20	S	S	S 333-20	S	S	S4	S4	S4

15:35

X100

SMARTX Sulfide Microplate Template

Name: Melanie

Plate ID: _____

Date: 7/28/26

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S	S	S	S	S	S	S	S	S
B	Std 1	Std 1	Std 1	S	S	S	S0	S0	S0	S	S	S
C	Std 2	Std 2	Std 2	S	S	S	Dup	Dup	Dup	S	S	S
D	Std 3	Std 3	Std 3	S	S	S	S3	S3	S3	Spike	Spike	Spike
E	Std 4	Std 4	Std 4	S	S	S	S	S	S	Spike	Spike	Spike
F	Std 5	Std 5	Std 5	S	S	S	S	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S	S	S	S	S	S	Dup	Dup	Dup
H	S	S	S	S	S	S	S	S	S	S4	S4	S4

~~10:06~~

10:06

SMARTX Sulfide Microplate Template

SMARTX 2026 July

Name: melanie

Plate ID: plate 6

Date: 7/22/26

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S ⁶³³⁻ ₂₀	S	S	S ⁴⁴¹⁻ ₁₀	S	S	S ⁴⁴²⁻ ₁₀	S	S
B	Std 1	Std 1	Std 1	S ⁶³³⁻ ₄₀	S	S	S0	S0	S0	S ⁴⁴²⁻ ₂₀	S	S
C	Std 2	Std 2	Std 2	S ⁶³³⁻ ₈₀	S	S	S ⁴⁴¹⁻ _{Dup 10}	Dup	Dup	S ⁴⁴²⁻ ₄₀	S	S
D	Std 3	Std 3	Std 3	S ⁶³³⁻ ₁₂₀	S	S	S3	S3	S3	S ⁶³³⁻ _{Spike 10}	Spike	Spike
E	Std 4	Std 4	Std 4	S ⁴⁴⁰⁻ ₁₀	S	S	S ⁴⁴¹⁻ ₂₀	S	S	S ⁴⁴¹⁻ _{Spike 20}	Spike	Spike
F	Std 5	Std 5	Std 5	S ⁴⁴⁰⁻ ₂₀	S	S	S ⁴⁴¹⁻ ₄₀	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S ⁴⁴⁰⁻ ₄₀	S	S	S ⁴⁴¹⁻ ₈₀	S	S	S ⁴⁴²⁻ _{Dup 40}	Dup	Dup
H	S ⁶³³⁻ ₁₀	S	S	S ⁴⁴⁰⁻ ₁₂₀	S	S	S ⁴⁴¹⁻ ₁₂₀	S	S	S4	S4	S4

X100

11:12

SMARTX Sulfide Microplate Template

SMARTX 2026 July

Name: Melanie

Plate ID: plate 7

Date: 7/22/26

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S 442-120	S	S	S 450-40	S	S	S 451-80	S	S
B	Std 1	Std 1	Std 1	S 443-10	S	S	S0	S0	S0	S 451-120	S	S
C	Std 2	Std 2	Std 2	S 443-20	S	S	S 450-40 Dup	Dup	Dup	S 452-10	S	S
D	Std 3	Std 3	Std 3	S 443-40	S	S	S3	S3	S3	S 442-80 Spike	Spike	Spike
E	Std 4	Std 4	Std 4	S 443-80	S	S	S 450-120	S	S	S 450-120 Spike	Spike	Spike
F	Std 5	Std 5	Std 5	S 443-120	S	S	S 451-10	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S 450-10	S	S	S 451-20	S	S	S 452-10 Dup	Dup	Dup
H	S 442-80	S	S	S 450-20	S	S	S 451-40	S	S	S4	S4	S4

11:49

X100

SMARTX Sulfide Microplate Template

SMARTX 2026 July

Name: Melanie

Plate ID: plate 8

Date: 7/22/26

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S ⁴⁵²⁻ ₄₀	S	S	S ⁴⁶⁰⁻ ₁₀	S	S	S ⁴⁶¹⁻ ₁₀	S	S
B	Std 1	Std 1	Std 1	S ⁴⁵²⁻ ₈₀	S	S	S0	S0	S0	S ⁴⁶¹⁻ ₂₀	S	S
C	Std 2	Std 2	Std 2	S ⁴⁵²⁻ ₁₂₀	S	S	Dup ⁴⁶⁰⁻ ₁₀	Dup	Dup	S ⁴⁶¹⁻ ₄₀	S	S
D	Std 3	Std 3	Std 3	S ⁴⁵³⁻ ₁₀	S	S	S3	S3	S3	S ⁴⁵²⁻ ₂₀ Spike	Spike	Spike
E	Std 4	Std 4	Std 4	S ⁴⁵³⁻ ₂₀	S	S	S ⁴⁶⁰⁻ ₂₀	S	S	S ⁴⁶⁰⁻ ₂₀ Spike	Spike	Spike
F	Std 5	Std 5	Std 5	S⁴⁵³⁻₄₀	S	S0	S⁴⁶⁰⁻₄₀	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S⁴⁵³⁻₈₀	S	S	S⁴⁶⁰⁻₈₀	S	S	Dup ⁴⁶¹⁻ ₄₀	Dup	Dup
H	S ⁴⁵²⁻ ₂₀	S	S	S ⁴⁵³⁻ ₁₂₀	S	S	S ⁴⁶⁰⁻ ₁₂₀	S	S	S4	S4	S4

12:58

X100

SMARTX Sulfide Microplate Template

SMARTX 2026 July

Name: Melanie

Plate ID: plate 9

Date: 7/22/26

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 0	Std 0	S 461-120	S	S	S 463-40	S	S	S 453-80 X100	S	S
B	Std 1	Std 1	Std 1	S 462-10	S	S	S0	S0	S0	S	S	S
C	Std 2	Std 2	Std 2	S 462-20	S	S	S 463-40 Dup	Dup	Dup	S	S	S
D	Std 3	Std 3	Std 3	S 462-40	S	S	S3	S3	S3	S 461-80 Spike	Spike	Spike
E	Std 4	Std 4	Std 4	S 462-80	S	S	S 463-80	S	S	S 463-80 Spike	Spike	Spike
F	Std 5	Std 5	Std 5	S 462-120	S	S	S 463-120	S	S	S0	S0	S0
G	MC: 10ppt S5	MC: 10ppt S5	MC: 10ppt S5	S 463-10	S	S	S 453-40 X1	S	S	Dup	Dup	Dup
H	S 461-80	S	S	S 463-20	S	S	S 453-40 X100	S	S	S4	S4	S4

13:47

X 100 unless otherwise stated