

SMARTX: Porewater Ammonium

May 2025 Samples

2026-08-31

Contents

0.1	Import Data & Clean	3
0.2	Assessing standard Curves	3
0.3	Dilution Corrections - ensure the latest dilution is kept	6
0.4	Performance Check	7
0.5	Analyze the Check Standards	7
0.6	Analyze Blanks	8
0.7	Analyze Duplicates	9
0.8	Spikes	10
0.9	Matrix Effects	11
0.10	Unit Converted Data Column Added (mg/L to uM)	12
0.11	Sample Flagging - Within range of standard curve	12
0.12	Create Log for QAQC	12
0.13	Merge samples with Metadata & check all samples are present	12
0.14	Visualize Data	13
0.15	Export Processed Data	18

##Run Information

```
cat("Run Information: ") #lets you know what section you're in
```

Run Information:

```
#set the run date & user name
run_date <- "08/12/2025"
sample_year <- "2025"
sample_month <- "May"
user <- "Victoria Ellis"

#identify the files you want to read in
#read in as a list to accommodate multiple runs in a month
NH3_files <- c("Raw Data/SMARTX_May_2025_Run1.csv",
              "Raw Data/SMARTX_May_2025_Run2.csv",
              "Raw Data/SMARTX_May_2025_Run3.csv")

# Define the file path for QAQC log file - NO Need to change just check year
file_path <- "Raw Data/SEAL_SMARTX_QAQC_Log_2025.csv"
final_path <- "Processed Data/GCReW_SMARTX_Porewater_NH4_202505.csv"

#record any notes about the run or anything other info here:
run_notes <- "PE Check contaminated, NH3 high, PE Check changed after this run
             Moved metadata chunk to setup because 'here' package was needed"

#Add SMARTX Prefix if necessary to samples
Add_Prefix <- "YES"

#Collection_Dates = "Raw Data/Sample_Collection_Dates_2024.csv"

cat(run_notes)
```

```
## PE Check contaminated, NH3 high, PE Check changed after this run
##           Moved metadata chunk to setup because 'here' package was needed
```

##Setup

##Read in metadata and create similar sample IDs for matching to samples


0.1 Import Data & Clean

0.2 Assessing standard Curves

```
#Pull out standards data
```

```
## Assess Standard Curves
```


```
## 'geom_smooth()' using formula = 'y ~ x'
```




Test	Run_Date	R2_Flag	R2	Slope
Ammonia 2	08/08/2025	R2_Pass	0.9996549	2.050170
Ammonia 2	08/11/2025	R2_Pass	0.9999057	2.009256
Ammonia 2	08/12/2025	R2_Pass	0.9984853	1.934847

```
#Plot standards data
```

```
## 'summarise()' has grouped output by 'Test', 'Intercept', 'Slope', 'R_Squared',  
## 'Run_Date'. You can override using the '.groups' argument.
```





0.3 Dilution Corrections - ensure the latest dilution is kept

```
## [1] "Reran Samples: c(\"442-40\", \"442-80\", \"443-40\", \"443-80\", \"443-120\", \"450-40\", \"450-80\", \"450-120\", \"450-160\", \"450-200\", \"450-240\", \"450-280\", \"450-320\", \"450-360\", \"450-400\", \"450-440\", \"450-480\", \"450-520\", \"450-560\", \"450-600\", \"450-640\", \"450-680\", \"450-720\", \"450-760\", \"450-800\", \"450-840\", \"450-880\", \"450-920\", \"450-960\", \"450-1000\")"
```

```
## [1] "Diluted Samples: c(\"442-40\", \"442-80\", \"443-40\", \"443-80\", \"443-120\", \"450-40\", \"450-80\", \"450-120\", \"450-160\", \"450-200\", \"450-240\", \"450-280\", \"450-320\", \"450-360\", \"450-400\", \"450-440\", \"450-480\", \"450-520\", \"450-560\", \"450-600\", \"450-640\", \"450-680\", \"450-720\", \"450-760\", \"450-800\", \"450-840\", \"450-880\", \"450-920\", \"450-960\", \"450-1000\")"
```


```
## [1] "No Naming Issues Detected"
```

0.4 Performance Check

Test	Run_Date	PE_Flag	PE_Conc	PE_Target_Conc
Ammonia 2	08/08/2025	Performance Check NOT Within 25% - REASSESS	1.500342	1.034
Ammonia 2	08/11/2025	Performance Check Within 25% - PROCEED	1.277440	1.034
Ammonia 2	08/12/2025	Performance Check NOT Within 25% - REASSESS	1.443127	1.034

0.5 Analyze the Check Standards


Test	Run_Date	RSV_Flag	RSV	RSV_Cutoff
Ammonia 2	08/08/2025	RSV WITHIN RANGE - PROCEED	0.0225172	0.25
Ammonia 2	08/11/2025	RSV WITHIN RANGE - PROCEED	0.0225172	0.25
Ammonia 2	08/12/2025	RSV WITHIN RANGE - PROCEED	0.0225172	0.25




Test	Run_Date	CHK_Flag
Ammonia 2	08/08/2025	>60% of Checks Pass - PROCEED
Ammonia 2	08/11/2025	>60% of Checks Pass - PROCEED
Ammonia 2	08/12/2025	>60% of Checks Pass - PROCEED

0.6 Analyze Blanks

Test	Run_Date	BLK_Pct_Flag	Mean_Blkc_Conc	Quantile_25
Ammonia 2	08/08/2025	>60% of Blanks Pass - PROCEED	0.0052279	1.271356
Ammonia 2	08/11/2025	>60% of Blanks Pass - PROCEED	0.0072246	1.654718
Ammonia 2	08/12/2025	>60% of Blanks Pass - PROCEED	0.0098191	1.394393




0.7 Analyze Duplicates




Test	Run_Date	Dup_Flags
Ammonia 2	08/08/2025	>60% of Dups Pass - PROCEED
Ammonia 2	08/11/2025	>60% of Dups Pass - PROCEED
Ammonia 2	08/12/2025	>60% of Dups Pass - PROCEED

0.8 Spikes



Test	Run_Date	Spike_Flags
Ammonia 2	08/08/2025	>60% of Spikes have a Diff <50% - PROCEED
Ammonia 2	08/11/2025	>60% of Spikes have a Diff <50% - PROCEED
Ammonia 2	08/12/2025	>60% of Spikes have a Diff <50% - PROCEED


0.9 Matrix Effects



Test	Run_Date	Matrix_Flags
Ammonia 2	08/08/2025	Matrix Has Diff <50% - PROCEED
Ammonia 2	08/11/2025	Matrix Has Diff <50% - PROCEED
Ammonia 2	08/12/2025	Matrix Has Diff <50% - PROCEED

0.10 Unit Converted Data Column Added (mg/L to uM)

0.11 Sample Flagging - Within range of standard curve



##Add QAQC Flags

0.12 Create Log for QAQC

##Write QAQC Log

0.13 Merge samples with Metadata & check all samples are present


Merge Samples with Metadata

All sample IDs are present in data


0.14 Visualize Data

```
## Visualize Data
```

C3: Porewater NH3




C4: Porewater NH3




```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.  
## i Please use 'linewidth' instead.  
## This warning is displayed once every 8 hours.  
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was  
## generated.
```


C3: Porewater NH3 by Replicate



C4: Porewater NH3 by Replicate



0.15 Export Processed Data

““

#end